Medfor
Plastic Containers for Laboratories

www.medfor.co.uk

EU Declaration of Conformity (Directive 98/79/EC)

MANUFACTURER		
Name of Company	**Address**	**SRN**
Medfor Products Ltd	Unit 2, Gresham Industrial Estate, Eastern Road, Aldershot, GU12 4YD, UK.	GB-MF-000022996

AUTHORISED REPRESENTATIVE			
Name of Company	**Address**	**SRN**	**Telephone / Email**
Emergo Europe	Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands	NL-AR-000000116	+31.70.345.8570 emergoeurope@ul.com

PRODUCT IDENTIFICATION			
Product Name	**Product Name**	**Basic UDI_DI**	**Product REFs**
IVD Containers with Lids	IVD Containers with Lids, supplied sterile. Intended for general in-vitro use within the laboratory or as specimen receptacles.	506056603IVDF8	See details in Table 1 attached.

Intended use
The device is a sterile, single-use plastic receptacle intended to be used by trained healthcare professionals in clinical or laboratory settings for the collection, storage, and transportation of human biological samples, including urine, faeces, sputum, tissue for subsequent in vitro diagnostic testing. The device is not intended to participate in the analytical or measurement process but is critical for maintaining sample integrity prior to testing. It is not intended for self-testing or near-patient testing.

Risk Class	**Conformity Assessment Route**
General	98/79/EC Annex III excluding section 6

Statements
The devices have no measuring function and are not reusable surgical instruments.

Relevant Standards
Medfor Products Ltd manufactures the devices according to the applicable requirements of the relevant harmonised standards listed in Table 2 attached.

Medfor Products Ltd declares, under our sole responsibility, that the devices listed on the attached Device Schedule conform to the relevant provisions of EC Council Directive 98/79/EC dated 27 October 1998 on in vitro diagnostic medical devices and all relevant amendments.
Conformity assessment has been carried out in accordance with Annex III, excluding Section 6.

Medfor Products Ltd agrees to develop, implement and maintain a formally recognised Quality Management System to ensure continued adequacy and efficacy.

COMPANY REPRESENTATIVE /
PERSON RESPONSIBLE FOR REGULATORY COMPLIANCE : Rob Stephens

TITLE: Commercial Director

SIGNATURE:

PLACE: Aldershot, UK

DATE: 23 July 2026

Medfor Products Ltd
www.medfor.co.uk
sales@medfor.co.uk
Tel: +44 (0) 1252 371181
Aldershot, GU12 4YD, UK
VAT GB 709 2511 49
COMPANY NO. 3787886

TABLE 1: DEVICE SCHEDULE

Sub-group

Non-wrapped containers

Sterility Type

Sterile R

REF	Description	
BUB0270D	2.7L white PP bucket with fitted white PP tamper-evident snap-on lid, sterile R	5060566036498
BUB0570D	5.7L white PP bucket tightly fitted with white PP tamper-evident snap-on lid, sterile R	5060566036504
OWA0250D	250ml white PP jar fitted with white PP screw cap (wadded), sterile R	5060566036894
OWA0250E	250ml white PP jar fitted with white PP screw cap (wadded), sterile R, plain label	5060566036900
OWA0350D	350ml white PP jar fitted with white PP screw cap (wadded), sterile R	5060566036931
OWA0350E	350ml white PP jar fitted with white PP screw cap (wadded), sterile R, plain label	5060566036948
OWA0500D	500ml white PP jar fitted with white PP screw cap (wadded), sterile R	5060566036955
OWA0500E	500ml white PP jar fitted with white PP screw cap (wadded), sterile R, plain label	5060566036962
OWA0750D	750ml white PP jar fitted with white PP screw cap (wadded), sterile R	5060566036979
OWA0750E	750ml white PP jar fitted with white PP screw cap (wadded), sterile R, plain label	5060566036986
OWA1000D	1000ml white PP jar fitted with white PP screw cap (wadded), sterile R	5060566036993
OWA1000E	1000ml white PP jar fitted with white PP screw cap (wadded), sterile R, plain label	5060566037006
OWA2000D	2000ml white PP jar fitted with white PP screw cap (wadded), sterile R	5060566037037
OWA2000E	2000ml white PP jar fitted with white PP screw cap (wadded), sterile R, plain label	5060566037044
PJB0300D	300ml clear PET jar fitted with white plastic screw cap (wood-pulp wad), sterile R	5060566037310
PJB0300E	300ml clear PET jar fitted with white plastic screw cap (wood-pulp wad), plain label, sterile R	5060566037327
PJB0300L	300ml clear PET jar fitted with white plastic screw cap (wood-pulp wad), plain label and sterile seal, sterile R	5060980290551
PJC0300D	300ml clear PET jar fitted with metal screw cap (flow-seal), sterile R	5060566037419
PJC0300E	300ml clear PET jar fitted with metal screw cap (flow-seal), plain label, sterile R	5060566037426
PJD0250D	250ml clear PET jar fitted with white wadded urea formaldehyde screw cap, sterile R	5060980291459
PJD0250E	250ml clear PET jar fitted with white plastic screw cap (wood-pulp wad), plain label, sterile R	5060980291473
PJE0250D	250ml clear PET jar fitted with metal screw cap (flow-seal), sterile R	5060980291510
PJE0250E	250ml clear PET jar fitted with metal screw cap (flow-seal), plain label, sterile R	5060980291374
PNA0150D	150ml Trans PP jar fitted with white HDPE screw cap, sterile R	5060566037495
PNA0200E	200ml trans PP jar fitted with HDPE screw cap, sterile R, plain label	5060566037525
PNA0250D	250ml transparent PP jar fitted with white HDPE screw cap, sterile R	5060566037556
PNA0250E	250ml Trans PP jar fitted with white HDPE screw cap, sterile R, plain label	5060566037563
PNA0350E	350ml Trans PP jar with white HDPE cap, sterile R, plain label	5060566037594
PNA0500D	500ml transparent PP jar fitted with white HDPE screw cap, sterile R	5060566037624
PNA0500E	500ml Trans PP jar fitted with white HDPE screw cap, sterile R, plain label	5060566037631
PNA1000D	1000ml Trans PP jar fitted with white HDPE screw cap, sterile R	5060566037662
PNC0025D	25ml transparent polypropylene jar loosely fitted with white HDPE tamper-evident screw cap, sterile R	5060566037754
PNC0050D	50ml transparent PP jar with loosely fitted white HDPE tamper-evident screw cap, sterile R	5060566037785
PNC0100D	100ml transparent PP jar loosely fitted with white HDPE tamper-evident screw cap, sterile R	5060566037815
PNE0025D	25ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060566038003
PNE0025E	25ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060566038010
PNE0050D	50ml trans PP jar with white fitted HDPE t/e screw cap, sterile R	5060566038027
PNE0100D	100ml Trans PP jar tightly fitted with white HDPE tamper-evident cap, sterile R	5060566038034
PNE0150D	150ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060566038041
PNE0150E	150ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060566038058
PNE0200D	200ml Trans PP jar tightly fitted with white HDPE tamper-evident cap, sterile R	5060566038065
PNE0200E	200ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060566038072
PNE0250D	250ml transparent PP jar tightly fitted with white HDPE tamper-evident cap, sterile R	5060566038089
PNE0250E	250ml Trans PP jar tightly fitted with white HDPE tamper-evident cap, sterile R, plain label	5060566038096
PNE0350D	350ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060566038102
PNE0350E	350ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060566038119
PNE0500D	500ml transparent polypropylene jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060566038126

PNE0500E	500ml Trans PP jar tightly fitted with white HDPE tamper-evident cap, sterile R, plain label	5060566038133
PNE1000D	1000ml transparent PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060566038140
SSA0060D	60ml clear PS jar fitted with flow-seal metal screw cap, sterile R	5060566038874
SSA0060E	60ml clear PS jar fitted with flow-seal metal screw cap, plain label, sterile R	5060566038881
SSA0060F	60ml clear polystyrene jar fitted with flow-seal metal screw cap, hospital label, sterile R	5060566038898
SSA0060K	60ml clear PS jar fitted with flow-seal metal screw cap, sterile seal, sterile R	5060566038935
SSA0100D	100ml clear PS jar fitted with flow-seal metal screw cap, sterile R	5060566038959
SSA0100E	100ml clear PS jar fitted with flow-seal metal screw cap, plain label, sterile R	5060566038966
SSA0100M	100ml clear PS jar fitted with flow-seal metal screw cap, plain label and sterile seal, sterile R	5060566039000
SSB0250D	250ml clear PS jar fitted with wadded metal screw cap, sterile R	5060980290858
SSB0250E	250ml clear PS jar fitted with wadded metal screw cap, plain label, sterile R	5060980290872
SSB0250F	250ml clear PS jar fitted with wadded metal screw cap, hospital label, sterile R	5060980290896
SSC0060D	60ml clear PS jar fitted with white HDPE screw cap, sterile R	5060566039079
SSC0060E	60ml clear PS jar fitted with white HDPE screw cap, plain label, sterile R	5060566039086
SSC0060F	60ml clear PS jar fitted with white HDPE screw cap, hospital label, sterile R	5060566039093
SSD0100D	100ml clear PS jar fitted with white PP screw cap (EPE wad), sterile R	5060566039147
SSD0100E	100ml clear PS jar fitted with white PP screw cap (EPE wad), plain label, sterile R	5060566039154
SSD0100M	100ml clear PS jar fitted with white PP screw cap (EPE wad), sterile seal, sterile R	5060566039192
SSE0150D	150ml clear PS jar fitted with flow-seal metal screw cap, sterile R	5060566039208
SSE0150E	150ml clear PS jar fitted with flow-seal metal screw cap, plain label, sterile R	5060566039215
TPA0035D	35 ml transparent PP jar with fitted white PE screw cap, sterile R	5060566039413
TPA0035E	35 ml transparent PP jar with fitted white PE screw cap, plain label, sterile R	5060566039420
TPA0055D	55ml transparent PP jar fitted with white PE screw cap, sterile R	5060566039451
TPA0055E	55ml transparent PP jar fitted with white PE screw cap, plain label, sterile R	5060566039468
TPA0055Q	55ml transparent PP jar fitted with white PE screw cap, 'Buffered Formalin' label, sterile R	5060566039475
TPA0125D	125ml transparent PP jar fitted with white PE screw cap, sterile R	5060566039482
TPA0125E	125ml transparent PP jar fitted with white PE screw cap, plain label, sterile R	5060566039499
TPA0180D	180 ml transparent PP jar fitted with white PE screw cap, sterile R	5060566039505
TPA0180E	180 ml transparent PP jar fitted with white PE screw cap, plain label, sterile R	5060566039512
UVA0030D	30ml clear PS universal tube fitted with HDPE screw cap, sterile R	5060566039567
UVA0030E	30ml clear PS universal tube fitted with HDPE screw cap, plain label, sterile R	5060566039574
UVA0030F	30ml clear PS universal tube fitted with HDPE screw cap, hospital label, sterile R	5060566039581
WPA0500D	500ml white PP jar fitted with white HDPE screw cap, sterile R	5060566039659
WPA0500E	500ml white PP jar fitted with white HDPE screw cap, plain label, sterile R	5060566039666
WPA0700D	700ml white PP jar fitted with white HDPE screw cap, sterile R	5060566039697
WPA0700E	700ml white PP jar fitted with white HDPE screw cap, plain label, sterile R	5060566039703
WPA1250D	1250ml white PP jar fitted with white HDPE screw cap, sterile R	5060566039734
WPA1250E	1250ml white PP jar fitted with white HDPE screw cap, plain label, sterile R	5060566039741
WPC0050D	50ml white PP jar loosely fitted with white HDPE tamper-evident screw cap, sterile R	5060566039833
WPC0100D	100ml white PP tamper-turn jar loosely fitted with white HDPE tamper-evident screw cap, sterile R	5060566039871
WPC0300D	300ml white PP jar loosely fitted with white HDPE tamper-evident screw cap, sterile R	5060566039918
WPC0500D	500ml white PP jar with loosely fitted white HDPE tamper-evident screw cap, sterile R	5060566039949
WPC1250E	1250ml white PP jar loosely fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060566039994
WPE0050D	50ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060980290063
WPE0050E	50ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060980290070
WPE0100D	100ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060980290087
WPE0100E	100ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060980290094
WPE0300D	300ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060980290100
WPE0300E	300ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060980290117
WPE0500D	500ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060980290124
WPE0500E	500ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060980290131
WPE0500R	500ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R (with separate labels)	5060980290148
WPE0700D	700ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060980290155
WPE0700E	700ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, plain label, sterile R	5060980290162
WPE1250D	1250ml white PP jar tightly fitted with white HDPE tamper-evident screw cap, sterile R	5060980290179

WPE1250E	1250ml white PP jar tightly fitted with HDPE tamper-evident screw cap, plain label, sterile R	5060980290186

Sub-group
Wrapped containers
Sterility Type
Sterile R

REF	Description	GTIN (unit)
OWA0250D-K	250ml white PP jar fitted with white PP screw cap (EPE wad) (1 PC/pack), double wrapped, sterile R	5060566035460
OWA0500D-K	500ml white PP jar fitted with white PP screw cap (EPE wad) (1 PC/pack), double wrapped, sterile R	5060566035477
OWA0750D-K	750ml white PP jar fitted with white PP screw cap (EPE wad) (1 PC/pack), double wrapped, sterile R	5060566035484
OWA1000D-K	1000ml white PP jar fitted with white PP screw cap (EPE wad) (1 PC/pack), double wrapped, sterile R	5060566035507
OWA2000D-K	2000ml white PP jar fitted with white PP screw cap (EPE wad) (1 PC/pack), double wrapped, sterile R	5060566035521
OWF0500E-K	500ml opaque white PP jar separate with white PP screw cap (1 PC/pack), plain label double wrapped, sterile R	5060566035613
S114-K	90mm unvented petri dishes double wrapped (1 PC/pack), sterile R	5060566035743
S240	25ml polypropylene tube fitted with yellow screw caps (1 PC/pack), labelled, double wrapped, sterile R	5060566035866
SSA0060D-K	60ml clear PS jar fitted with flow-seal metal screw cap (1 PC/pack), double wrapped, sterile R	5060566036108
SSA0060F-K	60ml clear PS jar fitted with flow-seal metal screw cap, hospital label (1 PC/pack), double wrapped, sterile R	5060566036115
SSA0100D-K	100ml clear PS jar fitted with flow-seal metal screw cap (1 PC/pack), double wrapped, sterile R	5060566036122
SSA0100F-K	100ml clear PS jar fitted with flow-seal metal screw cap, hospital label (1 PC/pack), double wrapped, sterile R	5060566036139
SSB0250D-K	250ml clear PS jar fitted with wadded metal screw cap (1 PC/pack), double wrapped, sterile R	5060980290919
SSB0250F-K	250ml clear PS jar fitted with wadded metal screw cap, hospital label (1 PC/pack), double wrapped, sterile R	5060980290933
SSE0150D-K	150ml clear PS jar fitted with flow-seal metal screw cap (1 PC/pack), double wrapped, sterile R	5060566036160
SSE0150F-K	150ml clear PS jar fitted with flow-seal metal screw cap, hospital label (1 PC/pack), double wrapped, sterile R	5060566036177
TPA0035D-K	35ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped, Sterile R	5060980290537
TPA0055D-K	55ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped, Sterile R	5060566036429
TPA0125D-K	125ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped, Sterile R	5060566036443
TPA0180D-K	180ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped, Sterile R	5060980290797
TPE0035D-K	35ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped (+ client label), Sterile R	5060566036184
TPE0055D-K	55ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped (+ client label), Sterile R	5060566036191
TPE0055D-L10	55ml transparent PP jar with fitted PE screw cap (10 inner-packs of 1 PC/pack), tripple wrapped, Sterile R	5060566036207
TPE0125D-K	125ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped (+ client label), Sterile R	5060566036214
TPE0125D-L5	125ml transparent PP jar with fitted PE screw cap (5 inner-packs of 1 PC/pack), tripple wrapped, Sterile R	5060566036221
TPE0180D-K	180ml transparent PP jar with fitted PE screw cap (1 PC/pack), double wrapped (+ client label), Sterile R	5060566036238
TPF0055D-L5	55 ml transparent PP jar with fitted PE screw cap (5 inner-packs of 2 PCS/pack), tripple wrapped, Sterile R	5060566036245
TPF0125D-L4	125 ml transparent PP jar with fitted PE screw cap (4 inner-packs of 2 PCS/pack), triple wrapped, Sterile R	5060566036252
UVB0030D-K	30ml clear PS universal tube with fitted HDPE screw cap (1 PC/pack), double wrapped, Sterile R	5060566036269
UVB0030F-K	30ml clear PS universal tube with fitted HDPE screw cap, hospital label (1 PC/pack), double wrapped, Sterile R	5060566036276
WPC1250D-K	1250ml white PP jar fitted with white HDPE tamper-evident screw cap (1 PC/pack), double wrapped, sterile R	5060566036375
WPE0300E-K	300ml white PP jar with separate white HDPE tamper-evident screw cap, plain label (1 PC/pack), double wrapped, sterile R	5060566036382

TABLE 2: RELEVANT STANDARDS	
Applicable to all Products above:	
Standard	**Description**
EN ISO 13485:2016	Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes
EN ISO 14971:2019	Medical Devices - Application of Risk Magament to Medical Devices
EN ISO 15223-1:2021	Medical Devices - Symbols to be used with Medical Device Labels, Labelling and Information to be supplied – Part 1: General Requirements
EN ISO 20417:2021	Information supplied by the manufacturer of medical devices
EN ISO 6717:2021	In vitro diagnostic medical devices. Single use receptacles for the collection of specimens, other than blood, from humans
Applicable to gamma sterilised products above:	
Standard	**Description**
EN ISO 11137-1:2015	Sterilisation of health care products-Radiation-Part 1: Requirements for development, validation and routine control of a sterilisation process for medical devices
EN ISO 11137-2:2015	Sterilisation of health care products - Radiation - Part 2: Establishing the sterilization dose